



## Clonazepam [OVA] (DAG-WZ1020)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Clonazepam is a benzodiazepine drug having anxiolytic, anticonvulsant, muscle relaxant, and hypnotic properties. Clonazepam has an unusually long elimination half-life of 18–50 hours, making it generally considered to be among the long-acting benzodiazepines. Clonazepam is a chlorinated derivative of nitrazepam and therefore a chloro-nitrobenzodiazepine.
<b>Specificity</b>	N/A
<b>Species</b>	N/A
<b>Conjugate</b>	OVA
<b>Applications</b>	IA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Geigy 32883; Lexin; Neurotol; Sirtal; Telesmin