



Rosiglitazone[BSA] (DAG-WZ1049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Rosiglitazone is an antidiabetic drug in the thiazolidinedione class of drugs. Like other thiazolidinediones, the mechanism of action of rosiglitazone is by activation of the intracellular receptor class of the peroxisome proliferator-activated receptors (PPARs), specifically PPAR γ . Rosiglitazone is a selective ligand of PPAR γ , and has no PPAR α -binding action. Apart from its effect on insulin resistance, it appears to have an anti-inflammatory effect: nuclear factor kappa-B (NF κ B) levels fall and inhibitor (I κ B) levels increase in patients on rosiglitazone.

Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Rosiglitazon; Rosiglitazona; Rosiglitazonum