



Tiamulin[BSA] (DAG-WZ1055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tiamulin is a pleuromutilin antibiotic drug that is used in veterinary medicine particularly for pigs and poultry. Strict MRL has been established due to the potential side effect in human.
Specificity	N/A
Species	N/A
Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	2-(diethylaminoethyl)thioacetoxymutilin; 81723 HFU; Dynamutilin
----------	---