



Methadone [OVA] (DAG-WZ39100)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Methadone is a synthetic opioid, used medically as an analgesic and a maintenance anti-addictive for use in patients on opioids. Although chemically unlike morphine or heroin, methadone also acts on the opioid receptors and thus produces many of the same effects.

Specificity N/A

Species N/A

Conjugate OVA

Applications IA

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg, 5 mg

Buffer 0.01M PBS, pH 7.2

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords dl-Methadone; Dolophine; Phenadone; Amidone