



Acetyl Salicylic Acid [OVA] (DAG-WZ3990O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Aspirin, also known as acetylsalicylic acid (ASA), is a medication used to treat pain, fever, or inflammation. Specific inflammatory conditions which aspirin is used to treat include Kawasaki disease, pericarditis, and rheumatic fever. Aspirin is one of the most widely used medications globally, with an estimated 40,000 tonnes (44,000 tons) (50 to 120 billion pills) consumed each year.

Specificity N/A

Species N/A

Conjugate OVA

Applications IA

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg, 5 mg

Buffer 0.01M PBS, pH 7.2

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Aspirin