



## Phenytoin [OVA] (DAG-WZ4036O)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |  |
|----------------------------|--|
| <b>Antigen Description</b> | Phenytoin is an anticonvulsant drug which can be useful in the treatment of epilepsy. Phenytoin blocks voltage-sensitive sodium channels in neurons. This action leads to a delay in neuronal electrical recovery. |
| <b>Specificity</b>         | N/A  |
| <b>Species</b>             | N/A  |
| <b>Conjugate</b>           | OVA  |
| <b>Applications</b>        | IA   |
| <b>Format</b>              | Liquid   |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.  |
| <b>Size</b>                | 1 mg, 5 mg   |
| <b>Buffer</b>              | 0.01M PBS, pH 7.2  |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.  |

### BACKGROUND

|                 |   |
|-----------------|---|
| <b>Keywords</b> | Dilantin; Diphenylhydantoin; Epanutin; Phenytek |
|-----------------|---|