



Theophylline [OVA] (DAG-WZ4380O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Theophylline is a methylxanthine drug used in therapy for respiratory diseases such as COPD and asthma under a variety of brand names.
Specificity	N/A
Species	N/A
Conjugate	OVA
Applications	IA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	1,3-Dimethylxanthine; Elixophyllin; Theophyllin; Theolair; Nuelin; Theophylline anhydrous; Respbid; Theocin
-----------------	---