



Digoxin [OVA] (DAG-WZ603520)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Digoxin is a beta blocker drug originally derived from the foxglove plant, <i>Digitalis lanata</i> . Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).
Specificity	N/A
Species	N/A
Conjugate	OVA
Applications	IA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Digitalis; Lanoxin; Digoxin