



Haloperidol [OVA] (DAG-WZ8647O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Haloperidol is a butyrophenone antipsychotic drug that has strong central antidopaminergic activity. It is used in the treatment of schizophrenia, acute psychotic states and delirium.
Specificity	N/A
Species	N/A
Conjugate	OVA
Applications	IA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg, 5 mg
Buffer	0.01M PBS, pH 7.2
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Aloperidin; Bioperidolo; Brotopon; Butyrophenone antipsychotic drug; Dozic; Einalon S; Eukystol; Haldol; Halosten; Keselan; Linton; Peluces; Serenace; Serenase; Sigaperidol
----------	--