



# Inactivated *Morganella morganii* Culture Fluid (DAG-ZL0739)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Product Overview** This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.

**Conjugate** N/A

**Applications** Control

**Stability** 12 months

**Format** Liquid

**Concentration** CFU/mL (Batch dependent - please inquire should you have specific requirements)

**Size** 10 mL

**Buffer** Culture media

**Preservative** None

**Storage** Store at -20°C

## BACKGROUND

**Keywords** *Morganella morganii*; *M. morganii*