



Inactivated *Pseudomonas fluorescens* Culture Fluid (DAG-ZL1026)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Conjugate	N/A
Applications	Control
Stability	12 months
Format	Liquid
Concentration	CFU/mL (Batch dependent - please inquire should you have specific requirements)
Size	10 mL
Buffer	Culture media
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Pseudomonas fluorescens</i> ; <i>P. fluorescens</i>
-----------------	--