



Recombinant Human ZnT8 [His] (DAG-ZNT8)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	≥90%
Conjugate	His
Applications	EIA, ICA, CLIA
Molecular Weight	60.2 KDa
Format	Liquid
Size	1 mg
Buffer	0.02M PBS, PH7.4
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction Zinc transporter 8 (ZNT8) is a protein that in humans is encoded by the SLC30A8 gene. ZNT8 is a zinc transporter related to insulin secretion in humans. Certain alleles of the SLC30A8 gene may increase the risk for developing type 2 diabetes, but a loss-of-function mutation appears to greatly reduce the risk of diabetes.

Keywords SLC30A8; zinc transporter 8; ZnT-8; ZNT8