



## Triadimefon [KLH] (DAG027S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	N/A
Species	N/A
Conjugate	KLH
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Triadimefon is also known as Bayleton. It is a systemic fungicide in the triazole family of chemicals. Triadimefon is a fungicide used in agriculture to control various fungal diseases. As a seed treatment, it is used on barley, corn, cotton, oats, rye, sorghum, and wheat.
Keywords	1-(4-Chlorophenoxy)-3,3-dimethyl-1-(1H-1,2,4-triazol-1-yl)butan-2-one;Triadimefon