



Tetrodotoxin [KLH] (DAG035K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Tetrodotoxin is an extremely potent poison (toxin) found mainly in the liver and sex organs (gonads) of some fish, such as puffer fish, globefish, and toadfish (order Tetraodontiformes) and in some amphibian, octopus, and shellfish species.
Keywords	Tetrodotoxin; TTX; Neurotoxin