



## Lomefloxacin [KLH] (DAG049S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

### BACKGROUND

Introduction	Lomefloxacin is a fluoroquinolone antibiotic used to treat bacterial infections including bronchitis and urinary tract infections.
Keywords	Lomefloxacin