



Picolinate [HRP] (DAG075S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	HRP
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	Picolinate is a pyridinemonocarboxylate resulting from the removal of a proton from the carboxy group of picolinic acid. It is a conjugate base of a picolinic acid.
Keywords	picolinate; pyridine-2-carboxylate; 2-pyridinecarboxylate; pyridineat; 2-picolinate; Picolinic acid anion