



## Nalidixic Acid [KLH] (DAG093S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

### BACKGROUND

Introduction	Nalidixic acid (tradenames Nevigramon, Neggram, Wintomylon and WIN 18,320) is the first of the synthetic quinolone antibiotics. In the technical sense, it is a naphthyridone, not a quinolone: its ring structure is a 1,8-naphthyridine nucleus that contains two nitrogen atoms, unlike quinoline, which has a single nitrogen atom.
Keywords	Nalidixic Acid