



## Nadifloxacin [HRP] (DAG095S)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	HRP
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

### BACKGROUND

Introduction	Nadifloxacin is a fluoroquinolone antibiotic that is used topically. It is effective against the Gram-positive bacteria <i>S. aureus</i> and <i>P. acnes</i> , as well as Gram-negative bacteria.
Keywords	Nadifloxacin; nadifloxacin; Nadixa; Acuatim; S-Nadifloxacin; jinofloxacin