



Androstenedione(3) [HRP] (DAG1022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Androstenedione(3), HRP conjugate
Antigen Description	Androstenediones are a pair of isomeric steroids which differ at the 5-position. They are distinguished as 5 α -androstenedione and 5 β -androstenedione. In humans, 5 α -androstenedione is a metabolic precursor of both testosterone and estrone.
Species	N/A
Conjugate	HRP
Format	Concentrate
Size	0.5 ml
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Androstenediones are a pair of isomeric steroids which differ at the 5-position. They are distinguished as 5 α -androstenedione and 5 β -androstenedione. In humans, 5 α -androstenedione is a metabolic precursor of both testosterone and estrone.
Keywords	Androstenedione(3)