



## Huma.htmA (DAG107)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Protein purified from Human thyroid tissue
<b>Species</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid. The antigen is supplied in 0.05M sodium phosphate buffer, pH 7.5, with 0.15M sodium chloride.
<b>Concentration</b>	2.16 mg/ml
<b>Preservative</b>	None
<b>Storage</b>	2-8°C short term, -20°C long term

### BACKGROUND

<b>Introduction</b>	Iodination and coupling of the immunogenic tyrosines in thyroglobulin to yield the thyroid hormones T(3) and T(4).
<b>Keywords</b>	TPO; thyroid peroxidase; MSA; TPX; TDH2A; thyroperoxidase; thyroid microsomal antigen;