



Doramectin [HRP] (DAG115S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	HRP
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	Doramectin is a macrocyclic lactone derivative of ivermectin, an anthelmintic that potentiates glutamate- and GABA-gated chloride channel opening in nematodes.
Keywords	Doramectine; Doramectinum; Doramectin