



17α-Estradiol(3) [HRP] (DAG1181)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	17α-estradiol is an estradiol that is estra-1,3,5(10)-triene substituted by hydroxy groups at positions 3 and 17 (the 17α stereoisomer). It is a 17α-hydroxy steroid, a 3-hydroxy steroid and an estradiol.
Species	N/A
Conjugate	HRP
Applications	ELISA
Format	Concentrate
Concentration	Please refer to the vial label for the specific concentration.
Size	500 µL
Buffer	Supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution.
Storage	Store at 2-8 °C for up to 3 months or -20 °C for long term storage
Ship	Wet ice