



Methandriol(3b) [HRP] (DAG1234)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Methandriol(3b), HRP conjugate
Antigen Description	Methandriol is an orally active synthetic derivative of testosterone. It is an anabolic androgen, used illegally in animal husbandry. Methandriol is often administered as derivatives, such as methandriol dipropionate and methandriol diacetate. Methandriol is mainly excreted in urine as sulphate tracer.
Species	N/A
Conjugate	HRP
Format	Concentrate
Size	0.5 ml
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Methandriol is an androgen.
Keywords	Methandriol