



Ronidazole [HRP] (DAG1285)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Ronidazole, HRP conjugate
Antigen Description	Ronidazol, dimetridazole, metronidazole, and ipronidazole are all nitroimidazole drugs that have been extensively used as growth promoters and therapeutic anti-bacterial and antiprotozoal agents in food-producing animals.
Species	N/A
Conjugate	HRP
Format	Concentrate
Size	0.5 ml
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Ronidazole is an antiprotozoal agent used in veterinary medicine for the treatment of histomoniasis and swine dysentery. It may also have use for the treatment of Tritrichomonas foetus infection in cats.
Keywords	Ronidazole