



Coumaphos [KLH] (DAG128S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	Coumaphos is a synthetic organic thiophosphate, organochlorine compound, and organophosphate acetylcholinesterase inhibitor that is used as a pesticide. It is characterized as a tan or slightly brown crystalline solid with a mild sulfurous odor, and exposure occurs by inhalation, ingestion, or contact.
Keywords	COUMAPHOS; Asuntol; Coumafes; Asunthol; Baymix; Resitox