



Stanozolol [HRP] (DAG1292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Stanozolol, HRP conjugate
Antigen Description	Stanozolol (C ₂₁ H ₃₂ N ₂ O) is a synthetic steroid that has anabolic and androgenic properties. Unlike most injectable anabolic steroids, Stanozolol is not esterified and is sold as an aqueous suspension.
Species	N/A
Conjugate	HRP
Format	Concentrate
Size	0.5 ml
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Stanozolol (C ₂₁ H ₃₂ N ₂ O) is a synthetic steroid that has anabolic and androgenic properties. Unlike most injectable anabolic steroids, Stanozolol is not esterified and is sold as an aqueous suspension.
Keywords	Stanozolol