



Toltrazuril [HRP] (DAG1319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|--|
| Product Overview | Toltrazuril, HRP conjugate |
| Antigen Description | Toltrazuril only need to use the drug for two days in unique disinsection way. Toltrazuril has no untoward reaction with high Security, so it can achieve the optimum efficiency and do not affect the crops' growth and feed availability. Toltrazuril can mix use with any drugs. Toltrazuril can remarkably reduce the death rate, recover quickly and grow fast. |
| Species | N/A |
| Conjugate | HRP |
| Format | Concentrate |
| Size | 0.5 ml |
| Preservative | None |
| Storage | 2-8°C short term, -20°C long term |

BACKGROUND

| | |
|--------------|--|
| Introduction | Toltrazuril is a coccidiostat, an antiprotozoal agent that acts upon Coccidia parasites. |
| Keywords | Toltrazuril |