



Clorprenaline [KLH] (DAG131S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	1-(2-chlorophenyl)-2-isopropylaminoethanol is a member of the class of monochlorobenzenes that is chlorobenzene which is substituted by a 1-hydroxy-2-[(propan-2-yl)amino]ethyl group at position 2. It is a member of monochlorobenzenes, a member of ethanolamines and a secondary amino compound.
Keywords	CLORPRENALINE; Clorprenalina; 1-(2-chlorophenyl)-2-(propan-2-ylamino)ethanol; Clorprenalinum; Isoprophenamin