



Trazadone [HRP] (DAG1323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Trazadone, HRP conjugate
Antigen Description	Trazadone is a serotonin uptake inhibitor that is used as an antidepressive agent. It has been shown to be effective in patients with major depressive disorders and other subsets of depressive disorders. It is generally more useful in depressive disorders associated with insomnia and anxiety. This drug does not aggravate psychotic symptoms in patients with schizophrenia or schizoaffective disorders.
Species	N/A
Conjugate	HRP
Format	Concentrate
Size	0.5 ml
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Trazadone is an antidepressant of the serotonin antagonist and reuptake inhibitor (SARI) class. It is a phenylpiperazine compound. Trazadone also has anxiolytic and hypnotic effects. Trazadone has considerably fewer prominent anticholinergic (dry mouth, constipation, tachycardia) and sexual side effects than most of the tricyclic antidepressants.
Keywords	Trazadone