



Cephalosporin [BSA] (DAG135S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	The cephalosporins are a class of β -lactam antibiotics originally derived from Acremonium, which was previously known as "Cephalosporium". Together with cephamycins they constitute a subgroup of β -lactam antibiotics called cephems.
Keywords	Cephalosporin; Cephalosporium