



# Recombinant HCV Nucleocapsid (a.a. 2-119) [HRP] (DAG1383)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Recombinant HCV Core HRP Labeled Antigen containing amino acids 2-119 was expressed in E. coli and purified by proprietary chromatographic technique.
<b>Antigen Description</b>	The hepatitis C virus (HCV) core protein represents the first 191 amino acids of the viral precursor polyprotein and is cotranslationally inserted into the membrane of the endoplasmic reticulum. Hepatitis C virus (HCV) core is a viral structural protein; it also participates in some cellular processes, including transcriptional regulation. However the mechanisms of core-mediated transcriptional regulation remain poorly understood. Hepatitis C virus (HCV) core protein is thought to contribute to HCV pathogenesis through its interaction with various signal transduction pathways. In addition, HCV core antigen is a recently developed marker of hepatitis C infection. The HCV core protein has been previously shown to circulate in the bloodstream of HCV-infected patients and inhibit host immunity through an interaction with gC1qR.
<b>Species</b>	HCV
<b>Purity</b>	> 95% pure as determined by 10% PAGE (coomassie staining).
<b>Conjugate</b>	HRP
<b>Applications</b>	HCV Core antigen is suitable for ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
<b>Size</b>	100 µg, 500 µg, 1 mg
<b>Buffer</b>	25mM Tris-HCl pH-8, 1mM EDTA, 1.5M urea and 50%glycerol.
<b>Preservative</b>	None
<b>Storage</b>	2-8°C short term, -20°C long term

# BACKGROUND

Introduction	Hepatitis C Virus is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10, 000 nucleotides and encodes a single polyprotein of about 3, 000 amino acids. The polyprotein is processed by host cell and viral proteases into three major structural proteins and several non structural proteins necessary for viral replication. Several different genotypes of HCV with slightly different genomic sequences have since been identified that correlate with differences in response to treatment with interferon alpha.
Keywords	HCcAg; Core protein p19; HCV core antigen; HCV core protein; Hepatitis C Virus core protein; HCV-1 Core Ag; Hepatitis C Virus Core Antige, genotype 6a; Flaviviridae; Hepacivirus