



Bromobutanol [KLH] (DAG140S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA, LF
Format	Liquid
Concentration	Please refer to the vial label for the specific concentration.
Size	1 mg
Buffer	Supplied in PBS
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid multiple freeze-thaw cycles.

BACKGROUND

Introduction	Among the β -agonist drugs illegally used as repartitioning agents in food producing animals, clenbuterol (CBL) and bromobutanol (BBL) are frequently recovered during routine control by national authorities.
Keywords	Brombuterol;Bromobutanol;AC1MJ3D4