



Recombinant HCV Nonstructural Protein 5 [HRP, Biotin] (DAG1441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant HCV NS5 Biotin Labeled protein containing the HCV NS5 immunodominant regions was expressed in E. coli and purified by proprietary chromatographic technique.
Species	HCV
Purity	> 95% pure as determined by 10% PAGE (coomassie staining).
Conjugate	HRP, Biotin
Applications	HCV-NS5 HRP antigen is suitable for ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.
Size	100 µg, 500 µg, 1 mg
Buffer	50mM Tris. pH-8 and 5mM EDTA.
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Hepatitis C Virus is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10, 000 nucleotides and encodes a single polyprotein of about 3, 000 amino acids. The polyprotein is processed by host cell and viral proteases into three major structural proteins and several non structural proteins necessary for viral replication. Several different genotypes of HCV with slightly different genomic sequences have since been identified that correlate with differences in response to treatment with interferon alpha.
---------------------	--

Keywords

HCV NS5; Hepatitis C virus non structural protein 5; Hepatitis C virus nonstructural protein 5;
Hepatitis C Virus NS5; HCV NS5 Genotype 2a; Hepatitis C Virus NS5 Genotype 2a
