



Sulindac [HSA] (DAG168S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sulindac, covalently linked to human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Sulindac is a sulfinylindene derivative prodrug with potential antineoplastic activity. Converted in vivo to an active metabolite, sulindac, a nonsteroidal anti-inflammatory drug (NSAID), blocks cyclic guanosine monophosphate-phosphodiesterase (cGMP-PDE), an enzyme that inhibits the normal apoptosis signal pathway; this inhibition permits the apoptotic signal pathway to proceed unopposed, resulting in apoptotic cell death. (NCI04)
---------------------	--

Keywords

drug-HSA conjugate, drug-conjugate, sulindac HSA conjugate
