



Bovine Fetuin (DAG175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Native BovineFetuin is a glycoprotein of MW48.4kD, and purified from fetal calf serum.
Species	Bovine
Conjugate	Unconjugated
Applications	ELISA
Size	1 g
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.
Keywords	A gamma globin; G gamma globin; Gamma 1 globin; Gamma 2 globin; Gamma A hemoglobin; Gamma globin; Hb F Agamma; Hb F Ggamma; hbfagamma; hbfaggamma; HBG 1; HBG 2; HBG1; HBG2; HBGA; HBGR; Hemoglobin gamma 1 chain; Hemoglobin gamma 2 chain; Hemoglobin gamma A;