



Recombinant HIV type 1 Pol Protein [GST] (DAG1851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	pol (HIV-1) full-length recombinant protein with GST tag expressed in Escherichia coli .
Species	HIV
Purity	Ion exchange column and HPLC reverse phase column
Conjugate	GST
Applications	SDS-PAGE
Format	Lyophilized
Buffer	Lyophilized from 20 mM PBS, 20 mM NaCl 0.5 M, pH 7.8 (1 mM DTT, 8 M urea)
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with distilled water, store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	HIV-1 is the most common and pathogenic strain of the virus. Scientists divide HIV-1 into a Group M (major) and two or more minor groups.
Keywords	Human Immunodeficiency Virus; HIV Pol; HIV; HIV-2 Pol; HIV-2; HIV type 2 Pol