



# Mouse Immunoglobulin Heavy Constant Gamma 3 (DAG189)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	MouseIgG3 with a kappa light chain.
Species	Mouse
Conjugate	Unconjugated
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

## BACKGROUND

Introduction	There are four IgG subclasses (IgG1, 2, 3, and 4) in humans, named in order of their abundance in serum (IgG1 being the most abundant).
Keywords	IgG3