



## Human IgM (DAG191)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Human
Conjugate	Unconjugated
Applications	ELISA
Preservative	None
Storage	2-8°C short term, -20°C long term

### BACKGROUND

Introduction	Immunoglobulin M, or IgM for short, is a basic antibody that is produced by B cells. IgM is by far the physically largest antibody in the human circulatory system. It is the first antibody to appear in response to initial exposure to antigen. The spleen is the major site of specific IgM production.
Keywords	IgM; Immunoglobulin M