



Biotinylated Recombinant SIV gp140 (mac239) (a.a. 20-692) [His] (DAG1976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	6xHis tagged HIV-1 gp140 (SIV/mac239) protein (a.a.20-692)
Species	SIV
Purity	≥95%
Conjugate	His
Applications	antibody ELISA, protein labeling, etc
Format	Liquid
Concentration	0.5 mg/ml
Size	50 µg
Buffer	PBS with 40% glycerol
Preservative	None
Storage	After reconstitution, stable for a month at 4°C. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing.

BACKGROUND

Introduction	The HIV-1 envelope (Env) glycoproteins play an essential role in the virus replication cycle by mediating the fusion between viral and cellular membranes during the entry process. The Env glycoproteins are synthesized as a polyprotein precursor, gp160, that is cleaved by cellular proteases to the mature surface glycoprotein gp120 and the transmembrane glycoprotein gp41.
---------------------	--

During virus assembly the gp120/gp41 complex is incorporated as heterotrimeric spikes into the lipid bilayer of nascent virions.

Keywords

Simian immunodeficiency virus; SIV gp140; SIV; SIV Glycoprotein 140
