



Recombinant HCV type 1b Core Antigen (a.a. 1-191) [His] (DAG1987)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Molecular Weight	22 kDa
Format	Liquid or Lyophilized powder
Size	100 µg, 1 mg
Buffer	PBS with 8M Urea
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords	Hepatitis C virus; HCV; Core protein
-----------------	--------------------------------------

GENE INFORMATION

UniProt ID	BAJ21643
-------------------	----------