



Recombinant HCV E2 Protein (a.a. 383-663) [His] (DAG2001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	6xHis tagged HCV E2 protein (a.a.383-663)
Species	HCV
Purity	≥ 95%
Conjugate	His
Applications	WB, etc
Format	Each vial contains 100 µg of purified protein in PBS containing 0.1% BSA and 25% glycerol.
Concentration	1 mg/ml
Size	100 µg, 1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Hepatitis C virus (HCV or sometimes HVC) is a small (55–65 nm in size), enveloped, positive-sense single-stranded RNA virus of the family Flaviviridae. Hepatitis C virus is the cause of hepatitis C in humans.
Keywords	HCV E2; Hepatitis C virus E2