



## Progesteron [BSA] (DAG208)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	Progesteron contains progesterone 3-CMO conjugated to Bovine Serum Albumin (BSA).
Species	N/A
Conjugate	BSA
Applications	ELISA
Size	25 mg
Preservative	None
Storage	Store at -20°C.

### BACKGROUND

Introduction	<p>Progesterone (pregn-4-ene-3,20-dione; abbreviated as P4) is an endogenous steroid hormone involved in the menstrual cycle, pregnancy, and embryogenesis of humans and other species. It belongs to a group of steroid hormones called the progestogens, and is the major progestogen in the body. Progesterone is also a crucial metabolic intermediate in the production other endogenous steroids, including the sex hormones and the corticosteroids, and plays an important role in brain function as a neurosteroid.</p>
Keywords	P4; Progesterone