



# Ofloxacin [HSA] (DAG208S)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Ofloxacin, covalently linked to human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

## BACKGROUND

Introduction	Ofloxacin is a fluoroquinolone antibacterial antibiotic. Ofloxacin binds to and inhibits bacterial topoisomerase II (DNA gyrase) and topoisomerase IV, enzymes involved in DNA replication and repair, resulting in cell death in sensitive bacterial species. (NCI04)
Keywords	drug-HSA conjugate, drug-conjugate, Ofloxacin HSA conjugate