



Recombinant HIV type 1 gp140 (SHIV SF162P3) (a.a. 30-668) [His] (DAG2139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal 6xHis tagged gp140 (SHIV SF162P3) (HIV-1) protein (a.a. 30-668)
Species	HIV
Purity	≥ 95%
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc
Format	Liquid
Concentration	1 mg/ml
Size	100 µg, 1 mg
Buffer	PBS
Preservative	None
Storage	Stable for a month at 4°C. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing.

BACKGROUND

Introduction	The HIV-1 envelope (Env) glycoproteins play an essential role in the virus replication cycle by mediating the fusion between viral and cellular membranes during the entry process. The Env glycoproteins are synthesized as a polypeptide precursor, gp160, that is cleaved by cellular proteases to the mature surface glycoprotein gp120 and the transmembrane glycoprotein gp41.
---------------------	--

During virus assembly the gp120/gp41 complex is incorporated as heterotrimeric spikes into the lipid bilayer of nascent virions.

Keywords

Human Immunodeficiency Virus; HIV gp140; HIV; HIV-1 gp140; HIV-1; HIV type 1 gp140
