



Recombinant HIV type 1 (AE/A244) gp70 [His] (DAG2146)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal 6xHis tagged gp70(MLV)-V1V2(HIV-1/CN54) fusion protein containing V1V2 loops from HIV-1(Clade AE/A244). A244 is the strain used in RV144(The Thai trial).
Species	HIV
Purity	≥ 95%
Conjugate	His
Applications	For probing V1V2 antibodies in immune sera from human clinical trial and animal studies
Format	Liquid
Concentration	1 mg/ml
Size	100 µg, 1 mg
Buffer	PBS
Preservative	None
Storage	Stable for a month at 4°C. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing.

BACKGROUND

Introduction	The HIV-1 envelope (Env) glycoproteins play an essential role in the virus replication cycle by mediating the fusion between viral and cellular membranes during the entry process. The Env glycoproteins are synthesized as a polyprotein precursor, gp160, that is cleaved by cellular
---------------------	--

proteases to the mature surface glycoprotein gp120 and the transmembrane glycoprotein gp41.

Keywords

Human Immunodeficiency Virus; HIV gp70; HIV; HIV-1 gp70; HIV-1; HIV type 1 gp70
