



Isoprenaline [HSA] (DAG228S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Isoprenaline, covalently linked to human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Isoprenaline is a secondary amino compound that is noradrenaline in which one of the hydrogens attached to the nitrogen is replaced by an isopropyl group. A sympathomimetic acting almost exclusively on beta-adrenergic receptors, it is used (mainly as the hydrochloride salt) as a bronchodilator and heart stimulant for the management of a variety of cardiac disorders. It has a role as a sympathomimetic agent, a beta-adrenergic agonist, a bronchodilator agent and a cardiotonic drug. It is a member of catechols, a secondary amino
--------------	---

compound and a secondary alcohol.

Keywords

drug-HSA conjugate, drug-conjugate, Isoprenaline HSA conjugate
