



# Recombinant CMV Pp150 (a.a. 495-691) [His] (DAG2390)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	C-terminal 6xHis tagged Pp150 (UL32) immunodominant region (a.a. 495-691)
<b>Species</b>	CMV
<b>Purity</b>	≥ 90%
<b>Conjugate</b>	His
<b>Applications</b>	WB standard, antibody ELISA, etc
<b>Format</b>	Each vial contains 100 µg of lyophilized protein in 500mM NaCl, 50mM Phosphat buffer (pH 7.4), 200mM imidazole, 8M urea.
<b>Concentration</b>	1 mg/ml
<b>Size</b>	100 µg, 1 mg
<b>Preservative</b>	None
<b>Storage</b>	2-8°C short term, -20°C long term

## BACKGROUND

<b>Introduction</b>	pp150 (UL32) is a basic phosphoprotein, one of the matrix components. It has been described as highly immunogenic, apparently more so than any other of the HCMV structural proteins.
<b>Keywords</b>	CMV p150; UL32; Cytomegalovirus pp150; Cytomegalovirus UL32; CMV HHV 5; Human herpesvirus 5; Protein UL150; UL150; Herpesviridae; Betaherpesvirinae; Cytomegalovirus