



Recombinant CMV Pp150 (a.a. 862-1048) [His] (DAG2391)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal 6xHis tagged Pp150 (UL32) immunodominant region (a.a. 862-1048)
Species	CMV
Purity	≥ 90%
Conjugate	His
Applications	WB standard, antibody ELISA, etc
Format	each vial contains 100 µg of lyophilized protein in 500mM NaCl, 50mM Phosphat buffer (pH 7.4), 200mM imidazole, 8M urea.
Concentration	1 mg/ml
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	pp150 (UL32) is a basic phosphoprotein, one of the matrix components. It has been described as highly immunogenic, apparently more so than any other of the HCMV structural proteins.
Keywords	CMV p150; UL32; Cytomegalovirus pp150; Cytomegalovirus UL32; CMV HHV 5; Human herpesvirus 5; Protein UL150; UL150; Herpesviridae; Betaherpesvirinae; Cytomegalovirus