



Recombinant HAV VP1 (a.a. 1-300) [His] (DAG2415)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal 6xHis tagged VP1 region of HAV capsid protein (a.a. 1-300)
Species	HAV
Purity	> 95%, based on SDS PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea.
Concentration	N/A
Size	100 µg, 1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Hepatitis A virus (HAV) is the sole member of the Hepatovirus genus within the family Picornaviridae. The capsid of HAV encloses a single-stranded RNA genome of about 7.5 kb which is translated into a single polyprotein. The virion proteins VP1 to VP4 and the nonstructural proteins are generated from the polyprotein by a cascade of proteolytic cleavages. Only one protease, viral protease 3C, has been implicated in the nine protein scissions. Processing of the capsid protein precursor region generates a unique intermediate, PX (VP1-2A), which accumulates in infected cells and is assumed to serve as precursor to VP1
---------------------	---

found in virions, although the details of this reaction have not been determined. Capsid proteins VP1, VP2, and VP3 form a closed capsid enclosing the viral positive strand RNA genome. VP1 is a major viral antigen.

Keywords

Group IV; Picornaviridae; Hepatovirus; Hepatitis A virus; Hepatitis A virus VP1 protein; HAV VP1; infectious hepatitis and epidemic virus
