



Diclofenac [HSA] (DAG244S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Diclofenac, covalently linked to human serum albumin
Species	N/A
Purity	chromatography, dialysis
Conjugate	HSA
Applications	e.g. allergy test, immunoassays etc.
Molecular Weight	app. 65 kDa
Format	Liquid
Concentration	1 mg/mL, conjugation ratio is app. 1:10 - 1:40.
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Diclofenac is a non-steroidal anti-inflammatory drug (NSAID) taken to reduce inflammation and as an analgesic reducing pain in certain conditions. The name is derived from its chemical name: 2-(2,6-dichloranilino) phenylacetic acid. In the United Kingdom, India, Brazil and the United States, it may be supplied as either the sodium or potassium salt, in China most often as the sodium salt, while in some other countries only as the potassium salt. Diclofenac is available as a generic drug in a number of formulations. Over-the-counter (OTC) use is
--------------	--

approved in some countries for minor aches and pains and fever associated with common infections.

Keywords

drug-HSA conjugate, drug-conjugate, Diclofenac HSA conjugate
